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Bibliographie

1. Aslan M, Yüksel P, Polat E, Cakan H, Ergin S, Öner YA, Zengin K, Arıkan S, Saribas S, Torun MM, Kocazeybek B. The diagnostic value of Western blot method in patients with cystic echinococcosis. *New Microbiol.* 2011; 34:173-177.
2. Bauomi IR, El-Amir AM, Fahmy AM, Zalat RS, Diab TM. Evaluation of purified 27.5 kDa protoscolex antigen-based ELISA for the detection of circulating antigens and antibodies in sheep and human hydatidosis. *J Helminthol.* 2015;89:577-583
3. Bhutani N, Kajal P. Hepatic echinococcosis: A review. *Ann Med Surg (Lond).* 2018;36:99-105
4. Deininger S, Wellingshausen N. Evaluation of a new combined Western and line blot assay (EUROLINE-WB) for diagnosis and species identification of Echinococcus infection in humans. *GMS Infect Dis.* 2019;7:Doc01.
5. Fathi S, Jalousian F, Hosseini SH, Parsa H, Kordafshari S. A Study of Cross-Reactivity Between Recombinant EPC1 Antigen of Echinococcus granulosus in Serum from Patients with Confirmed Cystic Echinococcosis Infection and Other Parasitic Infections. *Am J Trop Med Hyg.* 2016;94:1313-1317
6. Gao CH, Wang JY, Shi F, Steverding D, Wang X, Yang YT, Zhou XN. Field evaluation of an immunochromatographic test for diagnosis of cystic and alveolar echinococcosis. *Parasit Vectors.* 2018;11:311
7. Hadj Rabia S, Benmoussa F, Benzaid A, Baz A. Hydatidosis: Preparation and evaluation of radiolabeled antigens and antibodies. *Exp Parasitol.* 2018;187:67-74.
8. Hernández-González A, Sánchez-Ovejero C, Manzano-Román R, González Sánchez M, Delgado JM, Pardo-García T, Soriano-Gálvez F, Akhan O, Cretu CM, Vutova K, Tamarozzi F, Mariconti M, Brunetti E, Vola A, Fabiani M, Casulli A, Siles-Lucas M. Evaluation of the recombinant antigens B2t and 2B2t, compared with hydatid fluid,

- in IgG-ELISA and immunostrips for the diagnosis and follow up of CE patients. *PLoS Negl Trop Dis*. 2018;12:e0006741.
9. Khanbabaie S, Riazi M, Chang CH, Yunus MH, Noordin R. Lateral flow dipstick antigen assay for human cystic echinococcosis. *Acta Trop*. 2019;190:171-176.
 10. Lissandrin R, Tamarozzi F, Piccoli L, Tinelli C, De Silvestri A, Mariconti M, Meroni V, Genco F, Brunetti E. Factors Influencing the Serological Response in Hepatic Echinococcus granulosus Infection. *Am J Trop Med Hyg*. 2016;94:166-171
 11. Pagnozzi D, Biossa G, Addis MF, Mastrandrea S, Masala G, Uzzau S. An easy and efficient method for native and immunoreactive Echinococcus granulosus antigen 5 enrichment from hydatid cyst fluid. *PLoS One*. 2014;9:e104962.
 12. Manzano-Román R, Sánchez-Ovejero C, Hernández-González A, Casulli A, Siles-Lucas M. Serological Diagnosis and Follow-Up of Human Cystic Echinococcosis: A New Hope for the Future? *Biomed Res Int*. 2015;2015:428205
 13. Reiterová K, Auer H, Altintaş N, Yolasmaz A. Evaluation of purified antigen fraction in the immunodiagnosis of cystic echinococcosis. *Parasitol Res*. 2014;113:2861-2867.
 14. Sarink MJ, Koelewijn R, Slingerland BCGC, Tielens AGM, van Genderen PJJ, van Hellemond JJ. Performance of the commercially available SERION ELISA classic Echinococcus IgG test for the detection of cystic echinococcosis in clinical practice. *J Helminthol*. 2019;93:636-639.
 15. Sarkari B, Rezaei Z. Immunodiagnosis of human hydatid disease: Where do we stand? *World J Methodol*. 2015;5:185-195
 16. Savardashtaki A, Mostafavi-Pour Z, Arianfar F, Sarkari B. Comparison of the Utility of Recombinant B8/2 Subunit of the Antigen B, Native Antigen, and a Commercial ELISA Kit in the Diagnosis of Human Cystic Echinococcosis. *Iran Biomed J*. 2019;23:246-252
 17. Santivañez SJ, Arias P, Portocarrero M, Rodriguez S, Gonzalez AE, Gilman RH, Gavidia CM, Garcia HH. Serological diagnosis of lung cystic hydatid disease using the synthetic p176 peptide. *Clin Vaccine Immunol*. 2012;19:944-947.
 18. Sedaghat F, Sadjjadi SM, Hosseini SV, Kazemian S, Sarkari B. Evaluation of a simple Dot-ELISA in comparison with counter-current immunoelectrophoresis for diagnosis of human hydatidosis. *Clin Lab*. 2011;57:201-205
 19. Siles-Lucas M, Casulli A, Conraths FJ, Müller N. Laboratory Diagnosis of Echinococcus spp. in Human Patients and Infected Animals. *Adv Parasitol*. 2017;96:159-257
 20. Souki N, El Khattabi W, Bopaka RG, Aichane A, Afif H. Apport de la sérologie hydatique dans le diagnostic positif de l'hydatidose. *Revue des Maladies Respiratoires*. 2015;32:16
 21. Stojkovic M, Junghans T. Cystic and alveolar echinococcosis. *Handb Clin Neurol*. 2013;114:327-34
 22. Tamer GS, Dündar D, Uzuner H, Baydemir C. Evaluation of immunochromatographic test for the detection of antibodies against Echinococcus granulosus. *Med Sci Monit*. 2015;21:1219-22
 23. Velasco-Tirado V, Alonso-Sardón M, Lopez-Bernus A, Romero-Alegría Á, Burguillo FJ, Muro A, Carpio-Pérez A, Muñoz Bellido JL, Pardo-Lledias J, Cordero M, Belhassen-García M. Medical treatment of cystic echinococcosis: systematic review and meta-analysis. *BMC Infect Dis*. 2018;18:306.
 24. Wang JY, Gao CH, Steverding D, Wang X, Shi F, Yang YT. Differential diagnosis of cystic and alveolar echinococcosis using an immunochromatographic test based on the detection of specific antibodies. *Parasitol Res*. 2013;112:3627-33.
 25. Wuestenberg J, Gruener B, Oeztuerk S, Mason RA, Haenle MM, Graeter T, Akinli AS, Kern P, Kratzer W. Diagnostics in cystic echinococcosis: serology versus ultrasonography. *Turk J Gastroenterol*. 2014;25:398-404
 26. Zait H, Achir I, Guerchani MK, Hamrioui B. Profil épidémiologique de 290 cas d'échinococose kystique humaine diagnostiqués au CHU Mustapha d'Alger (2006 à 2011). *Pathol Biol (Paris)*. 2013;61:193-198.